

FCT Cutter

Code: E.226.xxx

Usage

To prepare samples for the FCT (Flat Crush Test).

Device description

The FCT Cutter is equipped with premium blades for cutting precise, circular FCT samples out of corrugated board. With its four fixation pins, the device prevents sliding of the sample while cutting. Samples with a thickness of 10 mm can be easily cut. A foam plate ensures safe operation without damaging the desk underneath. To cut a sample, the sample material is placed onto the foam plate, the Cutter is put on top and the sample is fixed with the four pins. Pushing the lever clockwise around once, the sample is cut.

Models

Code	Sample size	Cutting diameter
E.226.100	100 cm ²	112.84 mm
E.226.060	60 cm ²	79.79 mm
E.226.064	10 in ² (64 cm ²)	90.60 mm
E.226.032	5 in ² (32 cm ²)	64.10 mm

For more information about customized solutions do not hesitate to contact your personal distributor or PTE Austria directly!

